


Lab 4: Compare MAX30102 PPG sensor and custom analog PPG module (PSOC6)

Introduction

This lab helps the user to acquire the PPG signal using two different sensor setups: using the analog conditioning stages which uses one photodetector and three different wavelengths (green, red and IR) or using the commercial MAXREFDES117 PPG module that includes the MAX30102 sensor with only two wavelengths (red and IR).

The PPG EduKit platform is shown below. The module can be used with the CY8CPROTO-063-BLE board using the bridge adaptor provided. The analog PPG module includes one RGB LED, a photodetector, a transimpedance amplifier, a high pass filter, a low pass filter and an amplification stage. In this lab the PPG signal is acquired from the last amplification stage, thus the board has to be configured with the proper component values.



Objectives


After completing this lab, you will be able to:

- Understand how to set up a PSOC Creator project
- Use and modify the PPG EduKit software libraries
- Configure microcontroller peripherals in PSOC Creator (ADC, Timer, IRQ, UART)
- Understand multiwavelength measurement and control acquisition wavelength
- Read data from MAX30102 sensor
- Use the PPG EduKit GUI application and visualize multiple wavelength signals

Procedure


This lab is separated into steps that provide information on the detailed instructions that follow. Follow these detailed instructions to progress through the lab.

The lab includes 6 primary steps: create a project using PSOC Creator for CY8CPROTO-063-BLE board, integrate the software libraries in the newly created project (LED driver, digital potentiometer, MAX30102 sensor), configure the peripherals in order to read the PPG signal, configure UART communication and interrupts, send and receive data over UART, convigure MAX30102 sensor, and finally run the PPG EduKit GUI application to visualize data and control acquisition wavelength




1 Creation of the Project


- Open PSoC Creator
- Go to File \rightarrow New \rightarrow Project




 \bullet Select the target device \to PSOC6 \to <Launch Device Selector>









• Select the CYBLE-416045-2 Device



• Select "Empty schematic" and click "Next"



- In the "Select target IDE(s)" window, click "Next"
- Create a new Workspace and the project name is Compare_sensors



2 Integrate PPG EduKit software libraries

The next step is to integrate the software drivers that are provided in the PPG EduKit package. As seen in the application note, some external components need to be interfaced using SPI or I2C. Libraries are provided to speed up and to facilitate the development of PPG related applications, such that the user can focus on PPG signal acquisition or other types of algorithms.


In order to plot the PPG signal using the PPG EduKit GUI the user has to integrate the following libraries: TLC5925 (LED driver library), AD5273 (digital potentiometer library), custom PPG, serial frame (UART), MAX30102 module, milliseconds and utils library. Each library consists of a source file (.c) and a header file (.h) and can be found in PPG EduKit package.

Utils library

The library imports all the common libraries for all the PPG EduKit libraries and defines a common error handler.

```
23
                              81
                                 void HandleError(void)
24
    #include <stdint.h>
                              82 🖂 {
    #include <stdbool.h>
                              83 📥
                                       /* Disable all interrupts. */
    #include <stdlib.h>
26
                              84
                                        disable irq();
    #include <string.h>
27
                              85
    #include <math.h>
28
                              86
                                      /* Infinite loop. */
29
                              87
                                       while(lu) {}
                              88 | 1
```

• Import the source file (utils.c) in CM4 (Core1) \rightarrow Source Files




• Repeat the process for the header file (utils.h). Import the file in CM4 (Core1) → Header Files

AD5273 driver


AD5273 is a 64-position, one-time programmable (OTP) digital potentiometer that employs fuse link technology to achieve permanent program setting and can be configured through an I2C-compatible interface.

• Import the source file (AD5273.c) in CM4 (Core1) \rightarrow Source Files




• Repeat the process for the header file (AD5273.h). Import the file in CM4 (Core1) \rightarrow Header Files


• Go to the **TopDesign** schematic. In the component catalog (right panel at the right of the screen), write **I2C** and drag and drop the component into the schematic.




• Configure the I2C component by double click on it. Rename the component as I2C_BUS, set the mode as master and enable TX/RX FIFO buffers.



ullet Go to Design Wide Resources \to Pins and assign the I2C SCL and SDA pins as follows:



• Build the project to check if there is any error.




• Take a look at the application programming interface (AD5273.c) and try to understand driver functions.


TLC5925 driver

TLC5925 low-power 16-channel constant-current LED sink drivers to contain a 16-bit shift register and data latches, leading to converted serial input data into a parallel output format. The serial data is transferred into the device via **SDI** line at every rising edge of the **CLK** line. **LE** line latch the serial data in the shift register to the output latch, and the **OE** line enables the output drivers to sink current.


- Import the source file (TLC5925.c) in CM4 (Core1) \rightarrow Source Files
- Repeat the process for the header file (TLC5925.h). Import the file in CM4 (Core1) \rightarrow Header Files
- Go to the **TopDesign** schematic. In the component catalog (right panel at the right of the screen), write **Analog pin** and drag and drop the component into the schematic.




• Configure the Pin component by double click on it.




• Repeat the step for SDI, OE and LE pins.



 \bullet Go to Design Wide Resources \to Pins and assign the pins as follows:



• Build the project to check if there is any error.




• Take a look at the application programming interface (TLC5925.c) and try to understand driver functions.


Custom PPG driver

The Custom PPG library is intended to include the user defined algorithms such as HR, SpO2 or digital filters. In this lab, the library includes only the ADC reading routine using an timer triggered ISR.


• Go to the **TopDesign** schematic. In the component catalog (right panel at the right of the screen), write **ADC** and drag and drop the component into the schematic. Search for **Analog pin** and drag and drop the component into the schematic.




• Configure the ADC component by double click on it.




• Configure the ADC sample mode as single shot mode.




• Configure the Pin component by double click on it.




 \bullet Go to Design Wide Resources \to Pins and assign the $\mathbf{AMP_IN}$ pin as follows:




• Connect the analog pin to the ADC (press W to place a wire).




- Search for **Clock** and drag and drop the component into the schematic.
- Configure clock to 1 MHz



- Search for **Timer counter** and drag and drop the component into the schematic.
- Configure the counter.



- Search for **Interrupt** and drag and drop the component into the schematic. Rename the component as **adcRead_isr**.
- Interconnect the components




- Import the source file (custom_ppg.c) in CM4 (Core1) \rightarrow Source Files
- \bullet Repeat the process for the header file (custom_ppg.h). Import the file in CM4 (Core1) \to Header Files
- Build the project to check if there is any error.
- Take a look at the application programming interface (custom_ppg.c) and try to understand driver functions.

Serial Frame library

This library is a UART wrapper used to interface the PSOC with the PPG EduKit GUI. UART Communication stands for Universal asynchronous receiver-transmitter. It is a dedicated hardware device that performs asynchronous serial communication. It provides features for the configuration of data format and transmission speeds at different baud rates.

- Go to the **TopDesign** schematic. In the component catalog (right panel at the right of the screen), search for **UART** and drag and drop the component into the schematic.
- Configure the component.



 \bullet Go to Design Wide Resources \to Pins and assign the UART pin as follows:




- \bullet Import the source file (serial_frame.c) in CM4 (Core1) \rightarrow Source Files
- Repeat the process for the header file (serial_frame.h). Import the file in CM4 (Core1) \rightarrow Header Files
- Build the project to check if there is any error.
- Take a look at the application programming interface (serial_frame.c) and try to understand driver functions.


Milliseconds library

Before including the MAX30102 library, the milliseconds library has to be configured and included. As the Arduino function **millis**, we needed such a type of library to be used in PSOC. Porting code from Arduino to PSOC requires sometimes creating new libraries, since the microcontroller access layer is different for both platforms.

- Search for **Timer counter** and drag and drop the component into the schematic.
- Configure the counter. With a clock source of 1MHz, a 1000 period represents 1kHz (1ms).



- Search for **Interrupt** and drag and drop the component into the schematic. Rename the component as **isr_1ms**.
- Interconnect the components, use the same clock source (1MHz).



The isr_1ms increments a counter every millisecond. In such a way, the functionality of the millis function is achieved. For instance, the sensor can be pooled for a period of time using the newly milliseconds library.

MAX30102 library

The MAX30102 is an integrated pulse oximetry and heart-rate monitor module. It includes internal LEDs, photodetectors, optical elements, and low-noise electronics with ambient light rejection. The MAX30102 provides a complete system solution to ease the design-in process for mobile and wearable devices. The original library is provided by SparkFun as an Arduino library [4]. The library used in this project is a ported version of the

original library, since PSOC6 and Arduino platforms are not compatible. The driver interface includes multiple functions and that can be found in MAX30102.h file.

```
void MAX30102 NextSample(void);
void MAX30102 SetFIFOAverage(uint8 t numberOfSamples);
void MAX30102 EnableFIFORollover(void);
void MAX30102_SetLEDMode(uint8_t mode);
void MAX30102_SetADCRange(uint8_t adcRange);
void MAX30102_SetSampleRate(uint8_t sampleRate);
void MAX30102 SetPulseWidth(uint8 t pulseWidth);
void MAX30102 SetPulseAmplitudeRed(uint8 t amplitude);
void MAX30102 SetPulseAmplitudeIR(uint8 t amplitude);
void MAX30102_SetPulseAmplitudeGreen(uint8_t amplitude);
void MAX30102 SetPulseAmplitudeProximity(uint8 t amplitude);
void MAX30102 EnableSlot(uint8 t slotNumber, uint8 t device);
void MAX30102_ClearFIFO(void);
void MAX30102_Setup(uint8_t powerLevel, uint8_t sampleAverage, uint8_t ledMode,
void MAX30102 SoftReset (void);
bool MAX30102_ReadChannels(uint16_t *channelsBuffer, bool readFIF0);
bool MAX30102 SafeCheck(uint8 t maxTimeToCheck);
```

- Import the source file (MAX30102.c) in CM4 (Core1) \rightarrow Source Files
- Repeat the process for the header file (MAX30102.h). Import the file in CM4 (Core1) → Header Files

3 Create the main program

The main program should merge together all the libraries in such a way to reach the goal of measuring the data from both sensors and to send it over UART. Also, the user can select the LED wavelength. The main function can be found in main_c4.c.

• Go to main_c4.c and include the header file for each software driver.

```
20
                                               INCLUDE FILES
21
    * 1) system and project includes
    * 2) needed interfaces from external units
    ^{\star} 3) internal and external interfaces from this unit
23
25
26 - /* @brief Include PSOC generated files */
    #include "project.h"
28
29 □ /* @brief Include custom libraries for PPG EduKit */
   #include "utils.h"
30
   #include "AD5273.h"
31
32 #include "TLC5925.h"
   #include "custom_ppg.h'
33
   #include "serial_frame.h"
35 #include "MAX30102.h"
```

• Declare local macros. These macro directives are used for MAX30102 configuration, as parameters for MAX30102_Setup function.

```
41
                                          LOCAL MACROS
43
46 /* @brief MAX30102 Options: 0=Off to 255=50mA */
   #define CFG LED BRIGHTNESS
                               (60)
48 /* @brief MAX30102 Options: 1, 2, 4, 8, 16, */
49 #define CFG_SAMPLE_AVERAGE (4)
50 /* @brief MAX30102 Options: 1 = Red only, 2 = Red + IR, 3 = Red + IR + Green */
  #define CFG LED MODE
                                (2)
52 /* @brief MAX30102 Options: 50, 100, 200, 400, 800, 1000, 1600, 3200 */
   #define CFG SAMPLE RATE
                               (3200)
54 /* @brief MAX30102 Options: 69, 118, 215, 411 */
55 #define CFG_PULSE_WIDTH (411)
56 /* @brief MAX30102 Options: 2048, 4096, 8192, 16384 */
57
   #define CFG_ADC_RANGE
                               (16384)
58
```

• Declare the global variables. The serial frame command is represented by the **frameParams_t** structure type. The variable is passed to the **createSerialFrame** function in order to create the desired data frame. The function returns the address location of the data frame buffer to be sent, and **gpFrame** pointer is used to point to that memory location. Receiving data from UART is done using an interrupt request. Few variables have to be declared for this purpose. **Power**_levelvariable is used to store the current value for MAX 30102 sensor.

```
GLOBAL VARIABLES
75
76 L
78 - /* @brief Serial frame command */
   frameParams_t frameParam;
81 - /* @brief Pointer to a serial frame newly created*/
   uint8 t *gpFrame = NULL;
82
83
84 /* @brief Variables used for UART reading inside ISR */
   volatile uint8 t rx data[SERIAL FRAME LENGTH MAX] = {0};
   volatile uint32 t data rx index = 0;
86
87
   volatile bool frameReceived = FALSE;
88
   volatile uint8_t uartError;
90 1/2 / Obrief Variable used to set current for MAX30102 */
   uint8_t power_level;
```

• Before the main function declaration, **ISR_UART** function has to be declared. This function is later linked as the ISR handler for every RX interrupt request. Inside the function it is checked if the source of interrupt is due to "FIFO not empty". The interrupt flag is cleared and the data is read into **rx_data** buffer. The data is added to the buffer until the frame terminator is detected, then the buffer index is reset and the **frameReceived** flag is set as true.

```
101 - /*
102
                                               GLOBAL FUNCTIONS
103
104 void ISR UART (void)
105 □ {
106 占
         /* Check for "RX fifo not empty interrupt" */
         if((Uart Printf HW->INTR RX MASKED & SCB INTR RX MASKED NOT EMPTY Msk ) != 0)
107
108 🛱
109
             /* Clear UART "RX fifo not empty interrupt" */
110
             Uart Printf HW->INTR RX = Uart Printf HW->INTR RX & SCB INTR RX NOT EMPTY Msk;
111
             /* Get the character from terminal */
112
             rx_data[data_rx_index] = Cy_SCB_UART_Get(Uart_Printf_HW);
113
             /* Update data received index */
114
115
             if(rx data[data rx index] == FRAME TERMINATOR 2)
116
117
                 data rx index = 0;
118
                 frameReceived = TRUE;
119
120
             else
121
             {
122
                 frameReceived = FALSE;
123
                 data rx index++;
124
125
         }
126
         else
127
         £
128
             uartError = 1:
129
130 | }
```

• The main function starts with the initialization of adcRead_isr interrupt by setting the priority and the afferent interrupt vector. Then, the interrupts are enabled and the UART peripheral is initialized. Reading data from PC is done by configuring the UART RX FIFO interrupt and setting ISR_UART function as the handler.

The first component that is initialized is the AD5273 digital potentiometer, which is used for setting the current for the LED driver. The AD5273_Init function initializes the I2C peripheral, sets the data rate and the clock frequency. The LED current is set to 5 mA by calling the TLC5925_SetCurrent_mA function. The function translates the current value to a 6 bit value that represents the digital potentiometer resistance and writes the value over I2C. Having the external resistance configured, the LED driver can be configured to enable the infrared LED by calling TLC5925_enableIR. The PPG signal can be measured

right after the LED is turned on. Therefore, the ADC peripheral is started and the timer used to trigger the ADC conversion is initialized and started.

The MAX30102 sensor is initialized by calling the setup function with the desired configuration values. The parameters can take only few values and have been described above.

```
134 int main (void)
135 ⊟ {
136 🖨
         /* Assign ISR routines */
137
         CUSTOM PPG AssignISR AdcRead();
        /* Enable global interrupts. */
138
139
         enable irq();
140
         /* UART initialization status */
141
142
         cy en scb uart status t uart status ;
         /* Initialize UART operation. Config and Context structure is copied from Generated source.
143
         uart_status = Cy_SCB_UART_Init(Uart_Printf_HW, &Uart_Printf_config, &Uart_Printf_context);
144
         if (uart status != CY SCB UART SUCCESS)
145
146
147
             HandleError();
148
         Cy_SCB_UART_Enable(Uart_Printf_HW);
149
150
         /* Unmasking only the RX fifo not empty interrupt bit */
151
152
         Uart_Printf_HW->INTR_RX_MASK = SCB_INTR_RX_MASK_NOT_EMPTY_Msk;
153
         /* Interrupt Settings for UART */
         Cy_SysInt_Init(&Uart_Printf_SCB_IRQ_cfg, ISR_UART);
154
155
         /* Enable the interrupt */
156
         NVIC EnableIRQ(Uart Printf SCB IRQ cfg.intrSrc);
157
158
159
         /* Init digital potentiometer */
160
         AD5273_Init();
         /* Set 5 mA current for the LED driver */
162
         TLC5925 SetCurrent mA(5);
163
         /* Enable IR LED */
164
        TLC5925 enableIR();
         /* Start ADC and ADC Conversion */
165
166
         ADC_Start();
167
         /* Start timer used for ADC reading */
         CUSTOM_PPG_InitAndStartTimer_AdcRead();
168
169
         MAX30102_Setup(CFG_LED_BRIGHTNESS, CFG_SAMPLE_AVERAGE, CFG_LED_MODE,
                        CFG_SAMPLE_RATE, CFG_PULSE_WIDTH, CFG_ADC_RANGE);
170
```

• In the while loop, the **frameReceived** flag is checked every iteration. If a **SET_WAVELENGTH** frame is received, the value of the desired current value and wavelength are set. The current value for MAX30102 has to be mapped to a value between 5 and 255 (0.4 - 50 mA).

```
176
         /* Main infinite loop */
177
         while(1)
178
179
              if(frameReceived == TRUE)
180
181
                  frameReceived = FALSE:
                  if (rx_data[1] == SET_WAVELENGTH)
182
183
184
                      power_level = (255 / 45) * (rx_data[3] - 5);
185
                      activeChannel = rx_data[2];
186
187
                      if(rx data[3] > 40)
188
                          TLC5925 SetCurrent_mA(40);
189
190
                      1
191
                      else
192
193
                          TLC5925_SetCurrent_mA(rx_data[3]);
194
195
                      switch (activeChannel)
196
197
                         case IR CHANNEL:
198
                         TLC5925_enableIR();
199
                         MAX30102_SetPulseAmplitudeIR(power_level);
200
                         break;
201
                         case RED CHANNEL:
202
                         TLC5925 enableRed():
                         MAX30102 SetPulseAmplitudeRed(power_level);
203
204
                         break;
205
                         case GREEN CHANNEL:
206
                         TLC5925 enableGreen();
207
                          ^{\prime *} no green channel for MAX30102 (only for MAX30101) ^{*\prime}
208
                         MAX30102_SetPulseAmplitudeGreen(power_level);
209
                         break;
210
211
212
```

• The MAX30102 sensor should be read calling MAX30102_ReadChannels function in a busy loop. The FIFO_Buffer is filled until the head pointer reaches the maximum buffer size. A sliding window is used between two consecutive channel readings (the for loop is used to shift the buffers old values). Sending the entire buffer over UART leads to delay in displaying the signals, due to large frame sizes. Only the new acquired samples are sent over UART and plotted.

```
217
             for (uint8_t i = FIFO_NUMBER_OF_SAMPLES - FIFO_NUMBER_OF_OVERLAPPING_SAMPLES; i < FIFO_NUMBER_OF_SAMPLES; i++)
219
                 FIFO Buffer[((activeChannel * FIFO NUMBER OF SAMPLES) + (i - (FIFO NUMBER OF SAMPLES - FIFO NUMBER OF OVERLAPPING SAMPLES)))]
220
                                      = FIFO_Buffer[((activeChannel * FIFO_NUMBER_OF_SAMPLES) + i)];
222
223
             samplesTaken = FIFO NUMBER OF OVERLAPPING SAMPLES;
225
             /* Pool sensor until the number of samples is equal to FIFO NUMBER OF SAMPLES */
226
             while (FIFO NUMBER OF SAMPLES != samplesTaken)
228
                 MAX30102 ReadChannels(&FIFO Buffer[0], TRUE);
229
231
             /* Prepare serial frame structure */
232
             memset(&frameParam, 0x00, sizeof(frameParam));
             frameParam.frameType = CHANNEL_DATA;
234
             frameParam.sensor = MAX SENSOR;
235
             frameParam.params.wavelength = activeChannel;
             frameParam.tissueDetected = TRUE;
             gpFrame = createSerialFrame(&FIFO_Buffer[(activeChannel * FIFO_NUMBER_OF_SAMPLES) + FIFO_NUMBER_OF_OVERLAPPING_SAMPLES],
237
238
                                         (FIFO NUMBER OF SAMPLES - FIFO NUMBER OF OVERLAPPING SAMPLES) * 2, &frameParam);
239
             /* Sens serial frame */
240
             sendFrame (gpFrame);
243 - 1
245 □ /* [] END OF FILE */
```


- Depending on the MAX30102 sampling rate, the MAX30102_ReadChannels call can take some time to execute, since the function pools the sensor for new data over I2C. This is why the data acquired using the custom sensor is sent over UART using an interrupt request.
- In **custom_ppg.c** the interrupt handler for ADC reading is declared. The interrupt handler starts with clearing the interrupt flag. In the first if branch, the conversion is started and one of the buffers is filled with new data. It should be noticed that the I2C communication is interrupted before starting the ADC reading. It has been seen that this improves signal quality, since the I2C lines generates noise that is coupled into the analog signals.

```
111 void CUSTOM_PPG_InterruptHandler_AdcRead(void)
112 🖂 {
113
         /* Clear the terminal count interrupt */
114
         Cy_TCPWM_ClearInterrupt(Timer_AdcRead_HW, Timer_AdcRead_CNT_NUM, CY_TCPWM_INT_ON_TC);
115
116
117
         if((CUSTOM PPG BUFFER LENGTH != CUSTOM PPG bufferHead) && (TRUE == bBufferProcessed))
118
119
120
             if(I2C_BUS_MasterGetStatus() != CY_SCB_I2C_MASTER_BUSY)
121
122
                I2C BUS MasterAbortWrite();
                I2C BUS MasterAbortRead();
123
124
125
             ADC StartConvert();
127
             if(ADC IsEndConversion(CY SAR RETURN STATUS) != 0)
128
             {
129
                 switch (activeChannel)
130
131
                     case IR CHANNEL:
                     CUSTOM PPG bufferIR[CUSTOM PPG bufferHead] = ADC GetResult16(ADC CHANNEL 0 INV AMP);
132
133
                     break:
134
135
                     CUSTOM PPG bufferRed[CUSTOM PPG bufferHead] = ADC GetResult16(ADC CHANNEL 0 INV AMP);
136
137
                     case GREEN CHANNEL:
                     CUSTOM PPG bufferGreen[CUSTOM PPG bufferHead] = ADC GetResult16(ADC CHANNEL 0 INV AMP);
138
139
140
141
142
                 CUSTOM_PPG_bufferHead++;
143
144
         }
         else
145
```

• If the buffer head reaches the limit, the data is sent over the UART. As for MAX30102 buffer, an overlapping window is used between two different readings.


```
145
         else
146
147
             bBufferProcessed = FALSE;
             for (uintl6_t i = CUSTOM_PPG_BUFFER_LENGTH - CUSTOM_NUMBER_OF_OVERLAPPING_SAMPLES; i < CUSTOM_PPG_BUFFER_LENGTH; i++)
149
150
                 CUSTOM PPG bufferIR[(i - (CUSTOM PPG BUFFER_LENGTH - CUSTOM NUMBER_OF_OVERLAPPING SAMPLES))] = CUSTOM PPG bufferIR[i];
CUSTOM_PPG_bufferRed[(i - (CUSTOM_PPG_BUFFER_LENGTH - CUSTOM_NUMBER_OF_OVERLAPPING_SAMPLES))] = CUSTOM_PPG_bufferRed[i];
151
                 CUSTOM_PPG_bufferGreen[(i - (CUSTOM_PFG_BUFFER_LENGTH - CUSTOM_NUMBER_OF_OVERLAPPING_SAMPLES))] = CUSTOM_PFG_bufferGreen[i];
153
             ^{\prime } /* Prepare serial frame structure */
155
156
             memset((void *) &frameParamCustom, 0x00, sizeof(frameParamCustom));
157
             frameParamCustom.frameType = CHANNEL_DATA;
158
             frameParamCustom.sensor = CUSTOM SENSOR;
             frameParamCustom.params.wavelength = activeChannel;
160
             frameParamCustom.tissueDetected = TRUE:
             switch (activeChannel)
162
163
                case IR_CHANNEL:
164
165
                166
                                           (frameParams_t*) &frameParamCustom);
167
                break;
                case RED_CHANNEL:
                169
171
                                            (frameParams t*) &frameParamCustom);
172
                break;
173
                case GREEN CHANNEL:
                pFrame = createSerialFrame((uint16_t *) &CUSTOM_PPG_bufferGreen[CUSTOM_NUMBER_OF_OVERLAPPING_SAMPLES], (CUSTOM_PPG_BUFFER_LENGTH - CUSTOM_NUMBER_OF_OVERLAPPING_SAMPLES) * 2,
174
176
                                           (frameParams t*) &frameParamCustom);
178
             sendFrame((uint8_t *) pFrame);
bBufferProcessed = TRUE;
180
181
             CUSTOM_PPG_bufferHead = CUSTOM_NUMBER_OF_OVERLAPPING_SAMPLES;
```

• Now that both buffers are sent over UART, build the project and program the target.



4 Run PPG EduKit GUI application

The PPG EduKit GUI is a Python program used to interface the PPG EduKit platform with different types of applications. The user interface is created using QT Designer and stored in app2_gui.ui file. The main window loads the requested application in the list. The application has two source files: application2.py and Plotter2.py. The process_serial_frames.py file is used to interface the embedded application that communicates over UART with the PC. If the data frame structure is changed (by adding more functionalities), the change should be included both in the serial_frame.c file and in process_serial_frames.py script.



To run the application, one can use just a simple terminal, but the recommendation is to use Spyder IDE or PyCharm.

- Open the command line and go to the location of the PPG EduKit GUI application.
- Execute the command: **pip install -r requirements.txt**. If the pip package is not installed, please refer to [5] in order to set up your environment.
- Open the project in Spyder/PyCharm.
- Run the main.py file. In the application list interface, select the COMPARE PPG SOLUTIONS application.







UART PPG PLOTTING COMPARE PPG SOLUTIONS TBD

• Select the COM port and click **Connect**. Cover each sensor's surface and visualize the volumetric changes in arterial blood measured by two different sensors! Change the wavelength and the current to compare the optimal operating conditions for both modules. Also, one can change the sampling rate of both sensors in code and visualize the signal at different sampling rates.

```
void setup() {
  PPG_EK_Peripherals periphList = {
   .oledDisplay = ENABLE_PERIPHERAL,
    .neoPixel = DISABLE_PERIPHERAL,
   .tempSensor = DISABLE_PERIPHERAL,
   .ppgSensor = DISABLE_PERIPHERAL,
   .read_TIA = DISABLE_PERIPHERAL,
   .read_HPF = ENABLE_PERIPHERAL,
   .read_LPF = ENABLE_PERIPHERAL,
   .read_AMP = ENABLE_PERIPHERAL
  };
  delay(1000);
  PPG_Shield.begin(&periphList, 70);
  delay(1000);
  Adafruit SH110X &display reference = PPG Shield.getHandler OLED();
  display_reference.clearDisplay();
  PPG_Shield.enableLed(IR_CHANNEL, 5, true);
  curSwitch[AMP] = digitalRead(OLED BUTTON A);
  curSwitch[HPF] = digitalRead(OLED_BUTTON_B);
  curSwitch[LPF] = digitalRead(OLED_BUTTON_C);
  activeChannel = AMP;
 displayFilterText(AMP);
}
```

References

- [1] Cypress, "PSoC 6 MCU: CY8C63x6, CY8C63x7 Datasheet", https://www.cypress.com/file/385921/, Nov. 2020
- [2] Cypress, "PSoC 6 MCU Code Examples with PSoC Creator", https://www.cypress.com/documentation/code-examples/psoc-6-mcu-code-examples-psoc-creator, Mar. 2020
- [3] https://www.cypress.com/file/137441/download
- [4] https://github.com/sparkfun/SparkFun_MAX3010x_Sensor_Library
- $[5] \ \mathtt{https://phoenixnap.com/kb/install-pip-windows}$